



# M u l t i - P a n - C o a t e r

**GMPC**

GMPC III

GMPC II

GMPC I



**Glatt**<sup>®</sup>



# Glatt Multi-Pan-Coater

## GMPC

The new Glatt Multi-Pan-Coater (GMPC) is the newest development of the traditional Glatt Drum Coater. The GMPC is designed specifically for research and development for coating tablets with aqueous or organic solvent solutions.

As an option, it is also possible to coat cores with sugar syrup. The Multi-Pan has a compact construction and is suitable for GMP, features making it especially appropriate as a universal Coater.

All necessary components including the air intake and exhaust air system, are integrated into the machine housing. The Multi-Pan can be easily installed as a pilot machine.



The PLC control system records all necessary parameters, enabling the system to be scaled up to production plant size without any difficulty.

Due to the ease of the drum change, batch sizes from 2.0 to 56,0 litres can be processed.





# Addresses



## **Glatt Maschinen- & Apparatebau AG**

Kraftwerkstrasse 6  
4133 Pratteln / Switzerland  
Phone: +41 61 8 26 47 47  
Fax: +41 61 8 26 48 48  
eMail: info@glatt.com

## **Glatt customer service worldwide**



## **Glatt GmbH Process Technology**

Werner-Glatt-Straße 1  
79589 Binzen / Germany  
Phone: +49 7621 6 64 0  
Fax: +49 7621 6 47 23  
eMail: info@glatt.com

## **Glatt Ingenieurtechnik GmbH**

Nordstraße 12  
99427 Weimar / Germany  
Phone: +49 3643 47 0  
Fax: +49 3643 47 12 31  
eMail: info@glatt-weimar.de



## **Glatt Air Techniques Inc.**

20 Spear Road  
Ramsey, NJ 07446 / USA  
Phone: +1 201 8 25 87 00  
Fax: +1 201 8 25 03 89  
eMail: info@glattair.com

## **Glatt Norden ApS.**

Skøjtevej 27 - 31  
2770 Kastrup / Denmark  
Phone: +45 48 14 22 44  
Fax: +45 48 14 22 55  
eMail: info@glattnorden.dk

## **Glatt Pharmatech SASU**

Parc Technologique  
6, rue Louis Néel  
21000 Dijon / France  
Phone: +33 3 80 74 32 64  
Fax: +33 3 80 74 43 95  
eMail: info@glatt-pharmatech.fr

## **Glatt Protech Ltd.**

Swannington Road  
Cottage Lane Ind. Est.  
Broughton Astley  
Leicester LE9 6TU / Great Britain  
Phone: +44 1455 28 58 58  
Fax: +44 1455 28 55 10  
eMail: info@glatt-protech.co.uk

**www.glatt.com**

